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| STERILE MEDICAL DEVICE | DO NOT REUSE | SINGLE USE ONLY | LICENSE NO.: 25-440 |

- **Product Name:** RAY QUA
- **Model Name:** RAY QUA, RAY QUA HA25
- **Device Name:** Biomaterial for Tissue Restoration
- **Intended Use:** This product is intended for the temporary improvement of mild to moderate lateral canthal lines in adults through physical restoration by subcutaneous injection of sodium hyaluronate.

• Directions for Use

A. Preparations Before Use

1. Inspect the packaging for damage before use. Do not use the product if the packaging is damaged.
2. Check the expiration date indicated on the label. Do not use the product after the expiration date.
3. This product must only be used by a physician.
4. Allow the product to reach room temperature before injection.
5. Before treatment, inform the patient of the indications, contraindications, conditions for which the product is unsuitable, and potential adverse effects associated with the product.

B. Method of Use

1. Open the package containing the syringe.
2. Slowly remove the syringe cap.
3. Attach the needle securely to the syringe, and then remove the needle protective cap by pulling it horizontally.
4. Thoroughly disinfect the treatment area before injection.
5. Slowly inject the product into the skin. The amount injected may vary depending on the treatment area.
6. If the needle cannot be inserted smoothly, do not force it into the tissue. Withdraw the needle and perform the injection at another site.
7. After injection, gently massage the treatment area to ensure even distribution of the product.

C. Storage and Handling After Use

1. This product is a sterile, single-use medical device and must be discarded after use.
2. Do not resterilize the product.
3. Used syringes and needles must be disposed of as medical waste in accordance with applicable regulations.

D. Single-Use Device

This product is a sterile, single-use medical device. **Do not reuse.

• Precautions for Use

A. Warnings

1. Inspect the packaging for damage before use. Do not use the product if the packaging is damaged.
2. Check the expiration date indicated on the label. Do not use the product after the expiration date.
3. Do not reuse.
4. Do not resterilize.
5. Dispose of the used product in accordance with applicable regulations.

B. Precautions Regarding the Intended Patient Population, Including Age, Sex, and Health Status

1. Do not use this product in patients with known hypersensitivity to any component of the product, including sodium hyaluronate.
2. Do not use this product in areas previously treated with a non-absorbable filler or in areas that have recently undergone filler treatment.
3. This product is intended for the temporary improvement of mild to moderate lateral canthal lines in adults. Use in minors is prohibited.
4. Do not use this product in patients receiving non-steroidal anti-inflammatory drugs (NSAIDs), antiplatelet agents, anticoagulants, immunosuppressive agents, or similar medications.
5. Do not use this product in patients with allergic diseases, autoimmune diseases, sarcoid granulomatous disorders, or Osler’s endocarditis.
6. Use in pregnant or breastfeeding women is prohibited.
7. Do not use this product on damaged skin or in areas affected by inflammation or infection.
8. Do not administer excessive amounts of the product.

C. Precautions Regarding Adverse Events and Serious Complications That May Result from Use of the Medical Device or Improper Use

The physician must inform the patient of the potential medium- and long-term adverse effects associated with the use of this medical device.

1. Inflammatory reactions, including redness, swelling, and erythema, may occur at the injection site and may be accompanied by itching or pain, particularly when pressure is applied to the injection site. These reactions may persist for approximately one week.
2. Nodules or induration at the injection site, discoloration of the injection site, or an insufficient treatment effect may occur.
3. If inflammation persists for more than one week or if an adverse effect other than those described above occurs, the patient must be instructed to seek medical attention as soon as possible. The physician must provide appropriate treatment.
4. Any adverse effect other than those described above must be reported to the manufacturer or distributor.

D. General Precautions

1. Check the expiration date indicated on the label. Do not use the product after the expiration date.
2. Store away from direct sunlight.
3. Protect the product from impact or physical damage.

E. Interactions

Not applicable.

F. Use in Minors, Pregnant Women, and Breastfeeding Women

Use in minors, pregnant women, and breastfeeding women is prohibited.

G. Precautions for Application

1. Before use, inspect the packaging and confirm that the sterile barrier has not been compromised.
Do not use the product if the packaging or sterile barrier is damaged.
2. This product must be administered by a physician who has received sufficient training in the use of this product.
3. Before treatment, the physician must fully explain the indications, contraindications, and potential adverse effects of the product to the patient.
4. This product is intended for single use only. Do not reuse.
5. Avoid applying makeup for 12 hours after injection. Avoid prolonged exposure to sunlight, ultraviolet radiation, cold, or heat, including saunas and steam rooms, for two weeks after treatment.
6. The safety and effectiveness of long-term use beyond the period evaluated in clinical studies have not been established.

H. Precautions for the Prevention of Safety Incidents

1. Avoid use in vascular areas where vascular occlusion and subsequent tissue necrosis may occur.
2. Do not use the product until any infection or inflammation at the treatment site has been controlled or completely resolved.
3. The safety and effectiveness of this product for lip augmentation have not been established.
4. The safety of repeated injections of this product has not been established.
5. Inadvertent intravascular injection may result in serious adverse effects, including blindness. Use around the eyes, particularly in areas such as the glabella where the skin is thin and the risk of intravascular injection is high, is not recommended.
Exercise particular caution during injection.
6. Injection in patients with a history of herpetic eruptions may result in reactivation of the herpes infection.
7. The safety of this product has not been established in patients susceptible to keloid formation, hyperpigmentation, or hypertrophic scarring.

- **Shelf Life:** 36 months from the date of manufacture.
- **Storage Conditions:** Store at 2℃ to 25℃, Do not freeze. Protect from direct sunlight.
- **Packaging Unit:** 1 Syringe/box, 3 Syringes/box
- **Package Contents:** 1 syringe (2.5 mL) and 2 needles (30G) for the 1 syringe package; 3 syringes (2.5 mL each) and 6 needles (30G) for the 3 syringe package.
- **Manufacturer:** Re I Co., Ltd. 356-16 Gongdan-ro, Gunpo-si, Gyeonggi-do, Republic of Korea
- **Customer Service:** +82-70-4114-3802