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| STERILE MEDICAL DEVICE | DO NOT REUSE | SINGLE USE ONLY | LICENSE NO.: 25-159 |

Product Name: RAY NA

Model Name: RAY NA SHAPE

Device Name: Biomaterial for Tissue Restoration

Intended Use: This product is intended for the temporary improvement of moderate to severe nasolabial folds in adults through physical restoration by subcutaneous injection of cross-linked hyaluronic acid containing lidocaine.

Directions for Use

A. Preparations Before Use

1. Inspect the packaging for damage before use. Do not use the product if the packaging is damaged.
2. Check the expiration date indicated on the label. Do not use the product after the expiration date.
3. This product must only be used by a physician.
4. Allow the product to reach room temperature before injection.
5. Before treatment, inform the patient of the indications, contraindications, conditions for which the product is unsuitable, and potential adverse effects associated with the product.

B. Method of Use

1. Open the package containing the syringe.
2. Slowly remove the syringe cap.
3. Attach the needle securely to the syringe, and then remove the needle protective cap by pulling it horizontally.
4. Thoroughly disinfect the treatment area before injection.
5. Slowly inject the product into the skin. The amount injected may vary depending on the size and condition of the treatment area.
6. If the needle cannot be inserted smoothly, do not force it into the tissue. Withdraw the needle and perform the injection at another site.
7. After injection, gently massage the treatment area to ensure even distribution of the product.

C. Storage and Handling After Use

1. This product is a sterile, single-use medical device and must be discarded after use.
2. Do not resterilize or reuse the product.
3. Used syringes and needles must be disposed of as medical waste in accordance with applicable regulations.

D. Single-Use Device

This product is a sterile, single-use medical device. Do not reuse.

• Precautions for Use

A. Contraindications

Do not use this product in the following cases:

1. Patients with known hypersensitivity to any component of the product, including sodium hyaluronate.
2. Areas previously treated with a non-absorbable filler or areas that have recently undergone filler treatment.
3. Minors, as this product is intended for the temporary improvement of facial wrinkles in adults.
4. Patients receiving non-steroidal anti-inflammatory drugs (NSAIDs), antiplatelet agents, anticoagulants, immunosuppressive agents, or similar medications.
5. Patients with allergic diseases, autoimmune diseases, sarcoid granulomatous disorders, or Osler’s endocarditis.
6. Do not inject into blood vessels.
7. Do not use on damaged skin or areas affected by inflammation or infection.
8. Pregnant or breastfeeding women.
9. Use in minors is prohibited.
10. Do not administer excessive amounts of the product.

B. Warnings

1. Inspect the packaging for damage before use. Do not use the product if the packaging is damaged.
2. Check the expiration date indicated on the label. Do not use the product after the expiration date.
3. Do not reuse.
4. Do not resterilize.
5. Dispose of the used product in accordance with applicable regulations.

C. Adverse Effects

The physician must inform the patient of the potential medium- and long-term adverse effects associated with the use of this medical device.

1. Inflammatory reactions, including redness, swelling, and erythema, may occur at the injection site and may be accompanied by itching or pain, particularly when pressure is applied to the injection site. These reactions may persist for approximately one week.
2. Nodules or induration at the injection site, discoloration of the injection site, or insufficient treatment effect may occur.
3. If inflammation persists for more than one week or if an adverse effect other than those described above occurs, the patient must be instructed to seek medical attention as soon as possible. The physician must provide appropriate treatment.
4. Any adverse effect other than those described above must be reported to the manufacturer or distributor.

D. General Precautions

1. Avoid applying makeup for 12 hours after injection. Avoid prolonged exposure to sunlight, ultraviolet radiation, cold, or heat, including saunas and steam rooms, for two weeks after treatment.
 2. This product must be administered by a physician who has received sufficient training in the use of this product.
 3. Before treatment, the physician must fully explain the indications, contraindications, and potential adverse effects of the product to the patient.
 4. Before use, inspect the packaging and confirm that the sterile barrier has not been compromised.
- Do not use the product if the packaging or sterile barrier is damaged.
5. Check the expiration date indicated on the label. Do not use the product after the expiration date.
 6. Avoid use in areas where vascular occlusion and subsequent tissue necrosis may occur.
 7. Do not use the product until any infection or inflammation at the treatment site has been controlled or completely resolved.
 8. The safety and effectiveness of this product for lip augmentation have not been established.
 9. The safety of administering multiple injections of this product has not been established.
 10. Inadvertent intravascular injection may result in serious adverse effects, including blindness.
- Use around the eyes, particularly in areas such as the glabella where the skin is thin and the risk of intravascular injection is high, is not recommended. Exercise particular caution during injection.
11. Injection in patients with a history of herpetic eruptions may result in reactivation of the herpes infection.
 12. The safety of this product has not been established in patients susceptible to keloid formation, hyperpigmentation, or hypertrophic scarring.
 13. The safety and effectiveness of long-term use beyond the period evaluated in clinical studies have not been established.
 14. This product is intended for single use only. Do not reuse.
 15. Store away from direct sunlight.
 16. Protect the product from impact or physical damage.

• **Shelf Life:** 36 months from the date of manufacture.

• **Storage Conditions:** Store at 2°C to 25°C, Do not freeze. Protect from direct sunlight.

Packaging Unit: 1 Syringe/box

Package Contents: 1 syringe (1.0 mL) and 2 needles (25G)

Manufacturer: Re I Co., Ltd. 356-16 Gongdan-ro, Gunpo-si, Gyeonggi-do, Republic of Korea

Customer Service: +82-70-4114-3802